

Product Datasheet

ACTA2 antibody (HRP) (orb483041)

Description	ACTA2 antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Actin alpha
Target	alpha smooth muscle Actin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P62736
Dilution Range	WB=1:500-5000,IHC-P=1:100-500,IHC-F=1:500-1000,ICC=1:100
Expiration Date	12 months from date of receipt.