



Product Datasheet

CDK4 antibody (orb48321)


Description	Rabbit polyclonal antibody to CDK4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CDK4 (NP_000066.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	34kDa
Uniprot ID	P11802
Entrez	1019
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-100
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...