

Product Datasheet

Human CD25-PEPE antibody (PE) [Out of stock] (orb483573)

Description	Mouse monoclonal antibody against Human CD25-PEPE conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Target	IL2RA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	28 kDa
Uniprot ID	P01590
Dilution Range	Flow-Cyt=20ul/Test
Expiration Date	12 months from date of receipt.