


Product Datasheet

Integrin alpha 2 antibody (orb4839)

Description	Integrin alpha 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha 2 (1101-1181/1181 aa)
Target	Integrin alpha 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	127 kDa
Uniprot ID	P17301
Dilution Range	Flow-Cyt=1µg /test, WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of mouse kidney using Integrin alp...