


Product Datasheet

SIRP Alpha antibody (orb4840)

Description	SIRP Alpha antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD172a (151-250/504 aa)
Target	SIRP Alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P78324
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Western blot:Spleen(Mouse);Bones (Mouse)...