



## Product Datasheet

### MAP2K7 antibody (orb48447)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MAP2K7                                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | ICC, IF, WB                                                                    |
| <b>Immunogen</b>           | Recombinant fusion protein of human MEK7(NP_660186.1).                         |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                             |
| <b>Concentration</b>       | 1 mg/ml                                                                        |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                    |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | IgG                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                     |
| <b>MW</b>                  | 48KDa                                                                          |
| <b>Uniprot ID</b>          | <a href="#">O14733</a>                                                         |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200, ICC: 1:50-200, IF/ICC: 1:50-200, IP: 1:20-100 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |



Western blot analysis of extracts of var...



Immunofluorescence analysis of U2OS cell...