



Product Datasheet

DNA polymerase eta antibody (orb48495)


Description	Rabbit polyclonal antibody to DNA Polymerase theta
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human DNA polymerase eta(NP_006493.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80KDa
Uniprot ID	Q9Y253
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines.



IHC-P staining of human kidney tissue showing positive staining for DNA polymerase eta.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue showing positive staining for DNA polymerase eta.