



Product Datasheet

PRKCE antibody (orb48506)


Description	Rabbit polyclonal antibody to PRKCE
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PKC epsilon(NP_005391.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	Q02156
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200ICC1:50-200
Expiration Date	12 months from date of receipt.



IHC-P staining of human colon tissue using PRKCE antibody.



Immunohistochemical analysis of formalin-fixed rat spleen tissue using PRKCE antibody.



IHC-P image of rat spleen tissue using PRKCE antibody.