

Product Datasheet

isl1 antibody (orb48516)

Description	isl1 antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Insulin gene enhancer protein isl-1 protein (1-349AA)
Target	isl1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53405
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.