



Product Datasheet

S100A8 antibody (orb48527)


Description	Rabbit polyclonal antibody to S100A8
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human S100A8(NP_002955.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11KDa
Uniprot ID	P05109
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of MCF7 cell...



Western blot analysis of extracts of THP...