

Product Datasheet

GLP-1(1G9) antibody (PE) (orb485635)

Description	GLP-1(1G9) antibody (PE)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human GLP-1
Target	GLP-1(1G9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P01275
Expiration Date	12 months from date of receipt.