

Product Datasheet

IL1 beta antibody (PE) (orb485838)

Description	IL1 beta antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human IL-1beta
Target	IL-1 Beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01584
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, IF: 1:100-500, Flow-Cyt: 1µg/Test, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.