


Product Datasheet

Cdc25B (phospho-Ser186) antibody (orb4860)

Description	Cdc25B (phospho-Ser186) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Cdc25B around the phosphorylation site of Ser186 (AA(p-S)SS)
Target	phospho-Cdc25B (Ser186)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	P30306
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse brain tissue