


## Product Datasheet

### LEFTY1 antibody (orb48613)

|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| <b>Description</b>         | LEFTY1 antibody                                                         |
| <b>Species/Host</b>        | Rabbit                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                       |
| <b>Conjugation</b>         | Unconjugated                                                            |
| <b>Tested Applications</b> | ELISA, WB                                                               |
| <b>Immunogen</b>           | Recombinant Human Left-right determination factor 1 protein (223-260AA) |
| <b>Target</b>              | LEFTY1                                                                  |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol             |
| <b>Form/Appearance</b>     | Liquid                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.           |
| <b>Note</b>                | For research use only                                                   |
| <b>Isotype</b>             | IgG                                                                     |
| <b>Clonality</b>           | Polyclonal                                                              |
| <b>Uniprot ID</b>          | <b>075610</b>                                                           |
| <b>Dilution Range</b>      | WB:1:500-5000                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                         |



Western blot analysis of Rat spleen(lane...)