

Product Datasheet

GRIK3/GLR7 antibody (PE) (orb486171)

Description	GRIK3/GLR7 antibody (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human GRIK3/GLR7
Target	GRIK3
Preservatives	PBS (pH=7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	Q13003
Expiration Date	12 months from date of receipt.