

Product Datasheet

MGLUR1A/Metabotropic Glutamate Receptor 1a antibody (PE) (orb486174)

Description	MGLUR1A/Metabotropic Glutamate Receptor 1a antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human MGLUR1
Target	GRM1
Preservatives	PBS, Glycerol, rAlbumin
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13255
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.