

Product Datasheet

Glycine Receptor alpha 1 + alpha 2 antibody (PE)
[Out of stock] (orb486195)

Description

Rabbit polyclonal antibody

+ alpha 2 conjugated to PE

Species/Host

Rabbit

Conjugation

PE

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Glycine Receptor alpha 1 + alpha 2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

49kDa

Uniprot ID
Q64018, Q7TNC8, P22771, Q6DJV9, P23416, P23415, P07727
Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.