

Product Datasheet

DLL1 antibody (PE) (orb487191)

Description	Rabbit polyclonal antibody against DLL1 conjugated to PE
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human DLL1
Target	DLL1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76kDa
Uniprot ID	O00548
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.