

Product Datasheet

eIF4A3 antibody (PE) (orb487747)

Description	eIF4A3 antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human eIF4A3
Target	eIF4A3
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38919
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.