

Product Datasheet

TRPV2 antibody (PE) (orb488474)

Description	Rabbit polyclonal antibody against TRPV2 conjugated to PE
Species/Host	Rabbit
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human TRPV2
Target	TRPV2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86kDa
Uniprot ID	Q9WTR1
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.