

Product Datasheet

Intestinal alkaline phosphatase antibody (PE) [Out of stock] (orb489319)

Description	Intestinal alkaline phosphatase antibody (PE) [Out of stock]
Species/Host	Rabbit
Conjugation	PE
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ALPI
Target	ALPI
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09923
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, ICC: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.