



Product Datasheet

SIRT7 antibody (orb48950)


Description	Rabbit polyclonal antibody to SIRT7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SIRT7(NP_057622.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45KDa
Uniprot ID	Q9NRC8
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



IHC-P staining of human pancreas tissue ...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of HeLa cell...