



Product Datasheet

PIK3CB antibody (orb48953)


Description	Rabbit polyclonal to PIK3CB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Kinase p110 beta(NP_006210.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110KDa
Uniprot ID	P42338
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human colon tissue.



Western blot analysis of 293 cell and HepG2 cell lysates.