



Product Datasheet

CSRP1 antibody (orb48990)

Description	Rabbit polyclonal antibody to CSRP1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human CSRP1(NP_004069.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	P21291
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:20-50R, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of various cell li...



Immunoprecipitation analysis of 200ug ex...