



Product Datasheet

RNASE13 antibody (orb48992)


Description	Rabbit polyclonal to RNASE13
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human RNASE13(NP_001012264.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	Q5GAN3
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pr...



Immunohistochemical staining of human st...



Immunohistochemical staining of mouse lu...