

Product Datasheet

CHEK1 (phospho-Ser280) antibody (orb4903)

Description	CHEK1 (phospho-Ser280) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse CHEK1 around the phosphorylation site of Ser280 (AT(p-S)GG)
Target	phospho-CHEK1 (Ser280)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	014757
Expiration Date	12 months from date of receipt.