


Product Datasheet

CHEK1 (phospho-Ser317) antibody (orb4904)

Description	CHEK1 (phospho-Ser317) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CHEK1 around the phosphorylation site of Ser317 (SS(p-S)QP)
Target	phospho-CHEK1 (Ser317)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	014757
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...