

Product Datasheet

EGFR antibody (PE) (orb490695)

Description	EGFR antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human EGFR
Target	EGFR
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.