

Product Datasheet

Acetylated Lysine antibody (PE) [Out of stock] (orb490711)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody against Acetylated Lysine conjugated to PE |
| Species/Host | Mouse |
| Conjugation | PE |
| Tested Applications | ICC, IF |
| Immunogen | BSA conjugated Acetylated Lysine |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Lyophilized or Liquid |
| Concentration | - 100 µg (in 100 µl): 1 mg/ml |
| Storage | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Dilution Range | ICC: 1:50-200, IF: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |