

## Product Datasheet

**MKP1 (phospho-Ser296) antibody (PE) [Out of stock] (orb490726)**

**Description**

Rabbit polyclonal antibody against

to PE

**Species/Host**

Rabbit

**Conjugation**

PE

**Tested**

ICC, IF

**Applications**
**Immunogen**

KLH conjugated synthetic peptide derived from human Phospho-MKP1 (Ser296)

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

39kDa

**Uniprot ID**
**P28562, Q64623, P28563**
**Dilution Range**

ICC: 1:50-200, IF: 1:50-200

**Expiration Date**

12 months from date of receipt.