

Product Datasheet

CHK2 (phospho-ser516) antibody (orb4909)

Description	CHK2 (phospho-ser516) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CHEK2 around the phosphorylation site of Ser516 (QP(p-S)TS)
Target	Phospho-CHEK2 (Ser516)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	096017
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.