

Product Datasheet

Beta-Amyloid (1-40) antibody (PE) (orb491003)

Description	Mouse monoclonal antibody against Beta-Amyloid (1-40) conjugated to PE
Species/Host	Mouse
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human beta-Amyloid
Target	APP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	4.2kDa
Uniprot ID	P05067
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.