

Product Datasheet

MLKL (phospho-Ser358) antibody (PE) [Out of stock] (orb491019)

Description

Mouse monoclonal antibody against

to PE

Species/Host

Mouse

Conjugation

PE

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from human phospho-MLKL (Ser358)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

MW

54kDa

Uniprot ID

[Q9D2Y4](#), [Q8NB16](#)

Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.