

Product Datasheet

Tubulin, Alpha 4a antibody (PE) (orb491264)

Description	Tubulin, Alpha 4a antibody (PE)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Tubulin, Alpha 4a
Target	Tubulin, Alpha 4a
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68366
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.