

Product Datasheet

Myelin-oligodendrocyte glycoprotein antibody (PE) (orb491339)

Description

Rabbit polyclonal antibody against

conjugated to PE

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

PE

Tested Applications

ICC, IF

Immunogen

KLH conjugated synthetic peptide derived from human Myelin-oligodendrocyte glycoprotein

Target

MOG

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

24kDa

Uniprot ID
Q61885
Dilution Range

IF=1:100-500, ICC=1:100-500

Expiration Date

12 months from date of receipt.