



Product Datasheet

F11R antibody (orb49144)


Description	Rabbit polyclonal antibody to F11R
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human JAM-A/CD321/F11R(NP_058642.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37KDa
Uniprot ID	Q9Y624
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



IHC-P staining of human esophageal cancer...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of HeLa cell...