



Product Datasheet

PIK3R1 antibody (orb49152)

Description	Rabbit polyclonal antibody to PIK3R1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human Kinase p85 alpha(NP_852664.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85KDa
Uniprot ID	P27986
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of THP...



Immunohistochemical staining of rat brain...