

Product Datasheet

Human AIFM1 protein (orb49205)

Description	Recombinant human AIFM1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Form/Appearance	Supplied as a filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Glu121-Asp613
Source	E.coli
Uniprot ID	O95831
Expiration Date	6 months from date of receipt.