

Product Datasheet

HSV1 + HSV2 gD antibody (PE) (orb492164)

Description	Rabbit polyclonal antibody to Anti-HSV1 + HSV2 gD/PE
Species/Host	Rabbit
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human HSV1 + H
Target	HSV
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.