

Product Datasheet

Human BSSP4 protein (orb49224)

Description	Recombinant human BSSP4 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM HAc-NaAc, 150mM NaCl, 10% Glycerol, pH 4.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM HAc-NaAc, 150mM NaCl, 10% Glycerol, pH 4.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ala33-Ser317
Source	Human Cells
Uniprot ID	Q9GZN4
Expiration Date	6 months from date of receipt.