

Product Datasheet

Human BCAS2 protein (orb49225)

Description	Recombinant human BCAS2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 2mM DTT, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, 2mM DTT, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-T7&C-His
Note	For research use only
Application notes	T7 tag at the N-terminus, 6His tag at the C-terminus.
Protein Sequence	Ala2-Phe225
Source	E.coli
Uniprot ID	075934
Expiration Date	6 months from date of receipt.