

Product Datasheet

Human CA8 protein (orb49241)

Description	Recombinant human CA8 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 500mM NaCl, 1mM DTT, pH 8.5 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, 1mM DTT, pH 8.5 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ala2-Gln290
Source	E.coli
Uniprot ID	P35219
Expiration Date	6 months from date of receipt.