

Product Datasheet

LILRA3 antibody (PE) (orb492600)

Description	Rabbit polyclonal antibody against LILRA3 conjugated to PE
Species/Host	Rabbit
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human LILRA3/CD85e
Target	LILRA3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 2-8 °C for one year. Protect from light. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q8N6C8
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.