

Product Datasheet

Human CD40 Ligand protein (orb49280)

Description	Recombinant human CD40 Ligand protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Met113-Leu261
Source	E.coli
Uniprot ID	P29965
Expiration Date	6 months from date of receipt.