

Product Datasheet

Human DNA polymerase beta protein (orb49353)

Description	Recombinant human DNA polymerase beta protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 1mM EDTA, 50% Glycerol, pH 7.8.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 1mM EDTA, 50% Glycerol, pH 7.8.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ser2-Glu335
Source	E.coli
Uniprot ID	P06746
Expiration Date	6 months from date of receipt.