

Product Datasheet

Human STUB1 protein (orb49354)

Description	Recombinant human STUB1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Tyr303
Source	E.coli
Uniprot ID	Q9UNE7
Expiration Date	6 months from date of receipt.