

Product Datasheet

phosphor-PPP3CA (phospho-Ser469) antibody (PE)
[Out of stock] (orb493798)

Description

Rabbit polyclonal antibody

	(phospho-Ser469) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PPP3CA around the phosphorylation site of Ser469 (421-521/521aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	Q08209
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.