

Product Datasheet

Human FGF12 protein (orb49384)

Description	Recombinant human FGF12 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5mM EDTA, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 5mM EDTA, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Thr181
Source	E.coli
Uniprot ID	P61328
Expiration Date	6 months from date of receipt.