

Product Datasheet

Human FGF acidic protein (orb49387)

Description	Recombinant human FGF acidic protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 400mM NaCl, 1mM DTT, pH 8.0, The ED50 for this effect is less than 2 ng/ml, measured in a cell proliferation assay using BALB/c 3T3 cells.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM MOPS,100mMNa2SO4,1mM EDTA,pH7.9.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Phe16-Asp155
Source	E.coli
Uniprot ID	P05230
Expiration Date	6 months from date of receipt.