

Product Datasheet

Human FEN1 protein (orb49390)

Description	Recombinant human FEN1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50 mM Tris, 50mM NaCl, 1mM DDT, 10% Glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50 mM Tris, 50mM NaCl, 1mM DDT, 10% Glycerol, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Lys380
Source	E.coli
Uniprot ID	P39748
Expiration Date	6 months from date of receipt.